# ENDO-FLEX®

A MEDI-GLOBE COMPANY

## INSTRUCTION FOR USE

# ESOPHAGEAL MULTIBAND LIGATOR „MBL“ SU

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# MULTIBAND-LIGATUR-SET „MBL“ SU

**ENDO-FLEX GmbH**

Alte Hünxer Straße 115, D-46562 Voerde (Germany)

Phone: (0049) - (0)281 - 9 44 00 - 0

Fax: (0049) - (0)281 - 9 44 00 - 11

E-Mail: [sales@endo-flex.de](mailto:sales@endo-flex.de)

Website: [www.endo-flex.de](http://www.endo-flex.de)

Stand: 2019-12  
Current GA-0216

Page: 1 of 12  
Version: 5.0


*This document is only controlled in electronic format by the QM.*

**Attention:**

*This medical device may only be purchased by specialists, doctors and medical personnel and may only be used in accordance with these instructions for use and the defined area of application.*

**Description / Versions**

This device is composed of the following parts:

Luer-lock connection

Handle

Tracking shaft

Alignment arrows

Small Wheel

Luer-lock cap

Silicone sealing ring

Nylon tracking wire

Latex Ring

Penultimate ring in yellow

Band barrel

Introduction Catheter

**Scope**

This user manual is valid for the products listed below:

- MBL-01-06

**Important Note**

Please ensure you read this user manual carefully before each use and keep it somewhere that is easily accessible for users and/or appropriate specialist staff.

**Ensure you carefully read through the warnings. Patients, users or third parties may suffer serious injury as a result of the improper use of the products.**

## Content and Packaging

- 1 outer box
- 1 Multiband Ligator SU (Single Use) individually non-sterile packed
- 1 Instruction for use

## Patient Population

The patient population or target patient group is derived from the indication determined by the responsible doctor who diagnostically or therapeutically treats the patient as part of an endoscopic procedure (the leading procedure per se) according to the intended use of the medical device.

There are no known restrictions to the patient population or the target patient group.

### • Use of the product on minors

The user can use the product on minors in the case that the physiological and anatomical conditions of the patient permit the use of the product.

### • Use of the product on women who are pregnant or breastfeeding

The indication for the use of the product on women who are pregnant or breastfeeding must be narrowed down by the user based on the respective individual physiological and anatomical conditions of the patient.

## 1. Scope of Application

The products listed above must only be used by physicians with appropriate specialist training in gastroenterology. The products are intended solely for the medical sector as illustrated below and must therefore be used within an operating environment suitable for this purpose. It is essential that both the user as well as the appropriate specialist staff familiarize themselves with the instruments before the user makes use of them.

## 2. Application Period

These instruments are intended to be used uninterrupted for a period of maximum 60 minutes under regular conditions (MDD 93/42 EEC).

## 3. Intended Use

The Multiband Ligator is used via endoscope to ligate esophageal varices at and above the gastroesophageal junction.

### Indications

- Ligatur of esophageal varices

### Contraindication

The contraindications include without being limited to:

- cricopharyngeal or esophageal narrowing/stricture,
- tortuous esophageal diverticula,
- known or suspected esophageal perforation,
- asymptomatic rings or webs,
- coagulopathy.

Use of ligation bands is contraindicated in patients with a known hypersensitivity to latex.

## 4. Complications / Side Effects / Cross-reactions

Sedation during endoscopic examination increases the risk of hypoxemia, hypercapnia, hypotension, arrhythmias and aspiration due to reduced protective reflexes. Hypoxemia also occurs without sedation during endoscopic examinations due to the feed of the endoscope. Possible complications associated with the gastrointestinal endoscopy include and not limited to perforations, fever, aspirations, haemorrhage, infections, hypo-tension, allergic reactions to medication, cardiac arrhythmia or arrest, respiratory depression or arrest. Furthermore, complications that may occur with esophageal banding include and not limited to nausea, laryngeal, retrosternal pain, laceration, stricture formation, esophageal perforation, obstructions.

## 5. Warnings / Precautions

These instructions must be followed, as well as instructions from compatible components and hospital regulations for infection prevention, safe use, cleaning and sterilization. Failure to do so may result in serious injury to the patient and/or the user.

The following applies to the product:

- This device is delivered non-sterile and must NOT be sterilized before use!
- This device CANNOT be used with an over tube.
- For single use only! Do not reuse, reprocess or sterilize several times. Reuse, reprocessing or repeated sterilization of the instrument may affect its structural integrity and cause malfunction, resulting in contamination, infection and serious injury.

- If the instrument accidentally becomes dirty before treatment, it must be disposed of immediately! No cleaning agents may be applied.
- Do not use after the expiration date!
- All components should be carefully checked for compatibility and integrity before use. Do not use defective instruments! If defects occur, dispose of the instrument and replace it with a new one.
- Never use the product outside the recommended technical specifications (intended use).
- Never tamper with the structural conditions of the instrument, avoid kinks and other damage, and immediately discontinue use in the event of a malfunction!
- Wear protective clothing (gloves, goggles, gown, etc.) is absolutely necessary!
- Comparison of the technical data of the product with those of the endoscope used. The working channel diameter must be at least 0.2 mm larger than the outer diameter of the instrument.
- Never force instruments into the working channel!

It is vital the biopsy channel size of the endoscope is coordinated with compatible devices in order to obtain the best results during the procedure. The minimum biopsy channel diameter required is 2.8 mm and length of up to 180 cm. The outer diameter of the scope to be used must be from 8.5 mm up to 11.6 mm. Band ligation may not be effective when applied to small varices.

It is vital to ensure the understanding of the technical principles, clinical applications and associated risks with GI endoscopy as well as endoscopic banding before using this product. This product must only be used by or under direct supervision of physician thoroughly trained in therapeutic endoscopy as well as vessel banding.

Esophageal ligators are not intended for ligation of varices below the gastro-esophageal junction.

Banding should begin at the gastro-esophageal junction and proceed up the esophagus in an helicoidal manner, respecting a minimum of 2cm in between 2 ligatures. Passing the endoscope over an already ligated varices may dislodge the band from its position.

Routine endoscopic examination by the physician is recommended to confirm the diagnosis requiring treatment of esophageal varices or internal hemorrhoids prior to assembling the Multiband Ligator.

## 6. Compatibility

- Gastrosopes

## 7. Function Test

The medical devices must be checked with regard to the following aspects prior to use:

- Expiration date
- Undamaged packaging
- Damage on the product (cracks, bending, deformation)

Inspect products for immaculate surfaces, correct assembly and functionality.

## 8. Preparation / Application

- Check the instrument for reliable function and irregularities before the procedure.
- **If you notice any irregularities, replace the instrument with a new one!**
- When inserting into the working channel, the corresponding working channel and tube diameters must be observed.

### Procedure:


#### Mounting the MBL:

1. Unpack and check all parts of the system carefully. Ensure that all parts are in the package.
2. The nylon tracking wire of the Multiband Ligator system is already pre-mounted into an introduction catheter.
3. Straighten the endoscope. In case a biopsy valve is used, keep the lid open.
4. Insert the catheter into the biopsy channel until it comes to the distal tip of the scope and then release the catheter so that the Tracking wire is available for fixing it to the band barrel.
5. Fix the handle on the biopsy channel port.

**Hold the distal tip of the introduction catheter and remove the plastic part by slowly pulling it out of the endoscope until only tracking wire is left for the next fixing step.**

**Do not roll up the catheter on your hands to remove as it will block the inner tracking wire and then prevent the catheter to go out from the biopsy channel.**

6. Take the band barrel loop and connect it with the tracking wire loop properly as shown on image, making a flat knot:

**In order to tighten the knot, pull only the wires instead of pulling the whole cap as this could misfire a ring before the procedure.**

- Ensure that the handle has the alignment arrows facing each other.
- Using the small wheel, carefully and slowly pull the band barrel close to the tip of the scope.
- Mount the band barrel on the tip of the scope and ensure that the soft part is properly fixed on the scope in a way that it cannot fall off.
- Rotate the band barrel in order to prevent the tracking wire being in front of the endoscope camera.
- Slowly tension the tracking wire by turning the small wheel until you feel a small tension.

**Do not tension beyond this slight tension point to avoid shooting a ring.**

**If after inserting the handle of the device into the working channel entry port the handle is not stable enough or can move with soft movements of the scope, please remove the biopsy valve along with the handle of the device by 2cm up from the endoscope. Then press the handle further into the biopsy valve and place both handle and biopsy valve together into the endoscope working channel entry port again.**

#### Ligation of esophageal varices:

- Lubricate the scope and exterior part of the band barrel ensuring that no lubricant has entered the barrel.
- Roll back the big wheel 180° to permit that the tracking wires is not fully straightened so that no band is shot while inserting the endoscope into the patient.
- When in place, tension again the tracking wire using the small wheel to prepare for the shooting (be sure that the alignment arrows on the handle are facing each other).
- Locate the selected varix and aspirate it into the band barrel.
- To fire the band, slowly rotate 180° forward the big wheel. More than one ligation band for each varix might be needed to control acute bleeding.
- Release the suction button of the scope, insufflate air and then withdraw the scope slightly to release the ligated varice.
- Repeat the procedure until the last varix is banded.

**Note: The Luer-Lock connection is available to ensure irrigation of the scope biopsy channel if required. After removing the blue cap, attach the Luer-Lock connector to a syringe filled with sterile water and irrigate. After irrigation put back the blue cap to ensure good suction.**

**In case additional bands are needed, remove the scope and attach a new Multiband Ligator starting back from the mounting steps.**

#### Removing the MBL:

- Once the ligation procedure is completed, remove the scope from the patient.
- Dismantle the the Multiband Ligator as follows:
  - If any unused bands are still on the barrel shoot all remaining bands.
  - Remove the Multiband Ligator handle from the biopsy channel by pulling it, the remaining tracking wire attached with the handle should be removed also.
  - Remove the band barrel from the scope tip.
  - Dispose all parts per institutional guidelines for biohazardous medical waste.

**Note: If all bands have been used, you can remove the Multiband Ligator from the biopsy channel and continue using the scope and its biopsy channel while still inside the patient.**

## 9. Reprocessing and Sterilisation

These instruments are delivered in non-sterile condition and CANNOT EFFECTIVELY be cleaned, disinfected and sterilised after single use on account of the design which can no longer be removed and must be disposed of after single use.

## 10. Service Life of Products

The service life of the product is typically 2 years after the date of manufacture under normal conditions.

## 11. Liability and Warranty

As the manufacturer here, ENDO-FLEX shall not be liable for any damage or consequential damage arising from improper use or handling. This shall in particular apply to any use not in line with the defined purpose or failure to observe the preparation and sterilisation instructions and warnings. This shall also apply to repairs or modifications to the product made by individuals who are not authorised to do so by the manufacturer.

## 12. Service

Do not carry out any modifications to the product. If you have any concerns, complaints or comments regarding our products, please get in touch with us.

## 13. Return transport

Defective or non-conforming products must have gone through the entire reprocessing procedure prior to be sent back for repair/servicing. Please ensure you label the products accordingly with the note "hygienically safe" or "not decontaminated".

## 14. Transport and Storage Conditions

- Products may only be transported and stored in the packaging provided for this purpose.
- Products must be stored dry and protected from sunlight at room temperature.
- Do not place any objects on the storage packaging and the sterile barrier system!
- Do not store near aggressive chemicals.

## 15. Disposal

The products, packaging materials and accessories must all be disposed of in accordance with the nationally applicable regulations and laws. No specific instructions regarding this are given by the manufacturer.

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