
Kolondekompressionssonden SU

- DE
- EN
- FR
- IT
- ES
- NL
- CZ
- PT
- SE
- NO
- DK
- RU
- CN
- KR
- TU
- PL
- SA
- EL
- HR
- RO

DE	Kolondekompressionssonden SU (Einmalgebraucht)	Gebrauchsanweisung	3
EN	Colon Decompression Catheters (Single Use)	Instruction Manual	8
FR	Cathéter de décompression colon (Usage Unique)	Manuel d'utilisation	13
IT	Colon di decompressione del catetere (monouso)	Manuale di istruzioni	18
ES	Colón catéter de descompresión (De un solo uso)	Manual de instrucciones	23
NL	Colon decompressiekatheter (Eenmalig gebruik)	Handleiding	28
CZ	Dekompresní katetr pro tlusté střevo (Jednorázové použití)	Návod k použití	33
PT	Cateter de descompressão do cólon (Utilização única)	Manual de Instruções	38
SE	Colonic dekompression kateter (För engångsbruk)	Instruktionsbok	43
NO	Colonic dekompresjon kateter (Til engangsbruk)	Instruksjonshåndbok	48
DK	Kolonisk dekomprimeringskateter (Engangsbrug)	Betjeningsvejledning	53
RU	декомпрессионный катетер для толстой кишки (Одноразовая)	Руководство по эксплуатации	58
CN	结肠减压导管套件 (一次性使用)	说明书	63
KR	결장 압박 프로브 SU	사용 설명서	67
TU	Kolonik sıkıştırma probları su	Kullanım talimatları	72
PL	Sondy do kompresji okrężnicy SU	Instrukcje użytkowania	77
SA	مادخل سائل لتاميل عت (SU) ماعلا سل حمل		82
EL	Ανιχνευτές συμπίεσης παχέος εντέρου SU	Οδηγίες χρήσης	86
HR	Soniine kompresijske sonde SU	Uputa za uporabu	91
RO	Sonde de compresie colonică SU	Instrucțiuni de utilizare	96

Symbole

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Symbols

EN

	Medical device		Only for the attending physician
	Catalog number		Non-sterile
	Charge		Working length
	Attention! Read the instructions for use before use		Outer diameter
	Usable until		Guide wire maximum permissible diameter
	For single use only. Do not reuse		
	Do not sterilize again!		
	Keep dry		
	Protect from sunlight		
	Manufacturer		
	Date of manufacture		
	Do not use if packaging is damaged		

Copyright

The copyright to these operating instructions remains with the company ENDO-FLEX GmbH. All rights reserved, including those of translation. No part of this information may be reproduced in any form (print, photocopy, microfilm or any other process), either in whole or in part, without the written permission of the company ENDO-FLEX GmbH reproduced or processed, duplicated, translated or distributed using electronic systems. Violations may result in criminal penalties.

Limitation of liability

The ENDO-FLEX GmbH assumes no liability for personal injury, property damage, damage to the product or consequential damage resulting from failure to observe these operating instructions, from improper use of the product or from actions of unqualified personnel or personnel not trained by the ENDO-FLEX GmbH trained personnel. It is also prohibited to make unauthorized modifications or technical changes to the product.

Purpose of the instructions for use

These instructions for use are intended for medical personnel working with the product. It describes the structure, function and operation of the product. Read these instructions for use carefully and make sure you understand them and make them available to the healthcare professional

make!

Follow the instructions to ensure safe working. Also observe the information in the instructions of other devices used in combination with this product.

Do not follow instructions from unqualified third parties. The product may only be used by qualified medical personnel, who may be instructed by the ENDO-FLEX GmbH instructed medical personnel. The use of this product by healthcare professionals always involves weighing the individual patient's benefit against the patient's risk.

Always keep a copy of these instructions for use handy during use.

Classification

The Colonic compression probes SU is designated as a class I-Product according to Medical Device Regulation (MDR) - (2017/745) on medical devices.

Obligation to report

The ENDO-FLEX GmbH points out that the user is obliged to report all adverse incidents occurring during the use of this product and concerning it to the competent authority of his country and at the same time to ENDO-FLEX GmbH report.

Language

These and all other instructions for use are available at www.endo-flex.de as a PDF download.

1. Product description

A single lumen plastic catheter with radiopaque material. At the distal end over a length of 180 mm, there are six lateral openings to receive the air to direct it out of the colon. A funnel connector is attached to the proximal end. The Colonic Decompression Catheter is compatible with max. 0.52" guidewires.

These instructions for use are valid for the products listed below:

- 291140
- 291160
- 291180

1.1. Packing contents

- 1 outer carton
- 1 Colonic compression probes SU individually packaged
- 1 Instruction manual

1.2. Intended use

The Colonic compression probes SU are intended for the expulsion of air accumulation from the colon.

1.2.1 Application duration

The application time is specified for less than 60 minutes.

1.2.2 Lifetime

The service life of the product is 3 Years.

1.2.3 Indication

- Non-toxic megacolon
- Pseudo-obstructions
- Remaining air from previous colonoscopy

1.2.4 Contraindications

- Endoscopic decompression is contraindicated in the presence of bowel ischemia, bowel necrosis, and perforation. Colorectal cancer patients being treated or considered for treatment with anti-angiogenic therapy.

1.2.5 Potential complications/risks

- Under sedation during endoscopic examination, there is an increased risk of hypoxemia, hypercapnia, hypotension, arrhythmias, and aspiration due to decreased protective reflexes. Hypoxemia occurs even without sedation during endoscopic examinations due to the advancement of the endoscope.
- Possible complications associated with colonic decompression include perforation.

1.2.6 User target group

Trained specialists and professional staff.

1.3. Patient population

The patient or patient target group is derived from the indication of the responsible physician and the leading therapeutic procedure (interventional endoscopy for the removal of polyps and suspicious tissue in the gastrointestinal tract), who treats the patient according to the intended use of the medical device, diagnostically or therapeutically. Limitations to patient populations or patient target groups are not known.

1.3.1 Application for minors

The indication for the use of the product in minors must be narrowly defined by the user on the basis of the respective individual physiological and anatomical conditions.

1.3.2 Use in pregnant or breastfeeding women

The indication for use of the product in pregnant or lactating women must be narrowly defined by the user on the basis of the respective individual physiological and anatomical conditions.

1.4. Sterility

These single-use medical devices are supplied non-sterile. Sterilization is no longer possible and prohibited after use! The product may only be used once and must be disposed of after use in accordance with national and internal guidelines.

2. Safety instructions

The safety and warning information used in these operating instructions is based on ANSI Z535 standards.

DANGER	Warns of an imminent danger that will result in death or serious bodily injury.
WARNING	Warns of a potentially dangerous situation that could result in death or serious bodily injury.
CAUTION	Warns of a potentially dangerous situation that could result in moderate or minor bodily injury.
NOTE	Warns of a potentially dangerous situation that could lead to property and environmental damage.
SAFETY INSTRUCTION	Safety instructions (or equivalent signs) identify specific safety-related instructions or procedures.

2.1. General

NOTE

Never use the product outside the recommended technical specifications.

SAFETY INSTRUCTIONS

Protect the patient and yourself from contamination and injury by taking appropriate measures. Follow national and internal hygiene and safety guidelines. If in doubt, contact your local responsible person.

WARNING

Never tamper with the structural features of the instrument; abort the operation in the event of a malfunction.

WARNING

Do not use defective instruments! If defects occur, discard the instrument and replace it with a new one.

3. Application

3.1. Preparation

3.1.1 Work Preparation & Test

NOTE

Do not use products with expired expiration dates.

NOTE

If the packaging is already damaged before use, replace the instrument with a new one.

NOTE

If you notice any irregularities, replace the instrument with a new one.

WARNING

Using a defective instrument may cause damage to the endoscope and injury to patients or users.

3.1.2 Compatibility and accessories

- Guide wire up to max. .052".

3.2. Application

1. Insert the endoscope into the colon as far as possible or until you reach the caecum.
2. Check the position of the endoscope under X-ray or macroscopically.
3. Now guide a guide wire through the working channel of the endoscope until the distal tip emerges.
4. Now carefully remove the endoscope from the colon and make sure that the guide wire remains in the appropriate position. The position of the guide wire can be checked under X-ray vision.
5. Flush the colonic decompression catheter with sterile water or apply water-soluble lubricant.
6. Guide the colonic decompression catheter over the previously placed guide wire and check the correct position under X-ray vision.
7. Once the colonic decompression catheter is placed, carefully remove the guide wire.
8. Make sure that the funnel connector remains outside the body.

3.3. Preparation

WARNING

Reprocessing and sterilization is no longer possible and prohibited after use! The product may only be used once and must be disposed of after use in accordance with national and local guidelines.

4. Technical data

The list of all variants as well as the materials and individual parts used can be found in the Technical Product Data Sheet and is available on request from the ENDO-FLEX GmbH and/or on the website of ENDO-FLEX GmbH available.

5. Storage & Transport

This product may only be transported and stored in the packaging provided. Beyond that, no other special transport conditions are required.

Store the product at room temperature and protect it from direct sunlight and moisture. Do not store the instruments near chemicals, disinfectants or radioactive radiation.


SAFETY INSTRUCTIONS

Do not place any objects on the storage packaging and the sterile barrier system!

6. Disposal

After use, this product may present a biohazard. Disposal must be carried out in accordance with the regulations for clinically contaminated materials and national and regional laws and regulations.

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